

3D Accu<sup>i</sup>tomo 170




Thinking ahead.



Focused on life.

**THE PINNACLE OF 3D IMAGING:  
3D ACCUITOMO 170**



3D Accuitomo 170 is J. Morita's most advanced unit. Highly refined, it is the 4th generation of the Accuitomo product line. It produces the highest level of clarity that Morita offers with extremely low patient dosage.

#### **Reduced Dosage**

Dosage has been reduced 30% to 40% compared to the previous Accuitomo model.

#### **Four Imaging Modes**

A mode to serve every purpose. High Resolution Mode and High Fidelity Mode can be used for the highest quality of images. High Speed Mode reduces motion artifacts. Standard Mode can be used for both limited and broad fields of view.

#### **Nine Field of View (FOV) Sizes**

A wide FOV range, from 40 x 40 mm up to 170 x 120 mm, captures a couple teeth up to the entire head and neck area.

#### **Five Resolution Levels**

Select the voxel size that best suits your diagnostic needs: 80  $\mu\text{m}$ , 125  $\mu\text{m}$ , 160  $\mu\text{m}$ , 200  $\mu\text{m}$ , or 250  $\mu\text{m}$ .

#### **Zoom Reconstruction**

Use the original exposure data to zoom in on critical areas using voxel sizes as small as 80  $\mu\text{m}$ .

#### **DICOM Compatible**

Morita's i-Dixel 3D imaging software is DICOM 3.0 compliant.

#### **Viewing Software**

Allows viewing and manipulation of 3D data on computers without installation of 3D imaging software.

## INDUSTRY LEADING CLARITY WITH EXTREMELY LOW DOSAGE


**3D ACCUITOMO 170'S EXCEPTIONAL IMAGE QUALITY** can be attributed to several factors including the most advanced 3D sensors in the world and very high resolution at greater than 2 line pairs per millimeter.

### Super High Resolution, 80 $\mu\text{m}$ Images

Accuitomo's highly detailed images have a voxel size as small as 0.08 mm (or 80  $\mu\text{m}$ ). The slice width and pitch can be set between .08 and 3.2 mm. Depending on the diagnostic need, this unit offers five different resolution options to choose from: 80  $\mu\text{m}$ , 125  $\mu\text{m}$ , 160  $\mu\text{m}$ , 200  $\mu\text{m}$ , or 250  $\mu\text{m}$ .\*

### Extremely Low Patient Dose

3D Accuitomo 170 offers reduced dosage for increased patient protection. This unit's low dosage is rooted in how the machine acquires information. Through advanced engineering, the intensity of the X-rays has been optimized which decreases the overall level of emissions. Dosage has been reduced 30% to 40% compared to the previous Accuitomo model. Maximizing the efficiency of the machine has also resulted in sharper images of soft tissue with fewer artifacts. Now, even the largest field of view of the entire head and neck area is 25% less dosage than an FMX, F-speed film.\*\*



Even at the expanded  $\phi$  170 x H 120 mm size, 3D Accuitomo 170 maintains high resolution with low patient dose.

\* Depending on the size of the Field of View, some voxel sizes may not be possible.

\*\* Morita's effective dose calculated in accordance with ICRP 103 (ICRP 2007) FMX reference: Journal of the American Dental Association 2008; 139:1237-1243

## ZOOM RECONSTRUCTION

**THE ZOOM RECONSTRUCTION FEATURE CAN** generate a high resolution 80  $\mu\text{m}$  3D image for a region of interest from the larger size image already obtained.


Zoom Reconstruction is a feature unique to 3D Accuitomo. It always acquires data at 80  $\mu\text{m}$  the highest level of resolution. The user is able to choose a lower resolution for display, however, for large fields of view to reduce reconstruction time and file size. If a small region of interest is identified in the original scan, and more detail is needed in that area, Zoom Reconstruction will produce an 80  $\mu\text{m}$  image of that small region. This allows for higher level examination and diagnosis without retaking the image. No other unit offers this feature.

### Choice of Resolution Options


Voxel sizes of 80  $\mu\text{m}$ , 125  $\mu\text{m}$ , 160  $\mu\text{m}$  or 250  $\mu\text{m}$  can be selected as the resolution for each image area.

### Default Settings Include

- $\varnothing$  40 x H 40 mm: 80  $\mu\text{m}$
- $\varnothing$  60 x H 60 mm: 125  $\mu\text{m}$
- $\varnothing$  80 x H 80 mm: 160  $\mu\text{m}$
- $\varnothing$  100 x H 100 mm: 250  $\mu\text{m}$
- $\varnothing$  170 x H 120 mm: 250  $\mu\text{m}$



$\varnothing$  170 x H 120 mm, Voxel size: 250  $\mu\text{m}$



$\varnothing$  40 x H 40 mm, Voxel size: 80  $\mu\text{m}$

## FOUR IMAGING MODES

**MULTIPLE IMAGING MODES** allow selection of the scan time most appropriate for the clinical case, while keeping patient dose to a minimum.

### Standard Mode (Std)

360° scan: 17.5 seconds, 180° scan: 9.0 seconds.

Standard Mode offers high resolution images of exceptional clarity. This is the mode historically used in clinical case studies and research published on 3D Accuitomo units. Standard Mode is suitable for limited and wide views of temporal bone, paranasal sinus, maxilla and mandible, individual teeth, etc.

### High Resolution Mode (Hi-Res)

360° scan: 30.8 seconds, 180° scan: 15.8 seconds.

This is the highest resolution. Exposures are made at one-fourth the size of the detector pixels for the greatest spatial resolution. It is ideal for observation of delicate bone structures such as the ossicular chain. Available for 40 x 40 mm and 60 x 60 mm fields of view.

### High Fidelity Mode (Hi-Fi)


360° scan: 30.8 seconds, 180° scan: 15.8 seconds.

This mode has high density data to produce clear and sharp images. This mode is especially appropriate for performing zoom reconstructions.


### High Speed Mode (Hi-Speed)

360° scan: 10.5 seconds, 180° scan: 5.4 seconds.

High Speed mode reduces motion artifacts. It is a good choice for children or others who have difficulty remaining still. Available for 40 x 40 mm and 60 x 60 mm fields of view.




Std Mode  $\phi$  60 mm




Hi-Res Mode  $\phi$  60 mm







Hi-Fi Mode  $\phi$  170  $\times$  H 120 mm



Hi-Fi Mode  $\phi$  60 mm Zoom Reconstruction

# 3D IMAGE ACQUISITION




## How 3D Accuitomo 170 Acquires Images


The arm rotates 360° around the center of the exposure region for 18 seconds (Standard Mode) as the X-ray head emits a cone-shaped beam. The multiple projections created during the arm's rotation are converted to a digital signal by the flat panel detector and are transmitted to the computer. After any necessary supplemental or corrective processing, the digital information is converted into a three dimensional CT image using an image reconstruction algorithm, and a high resolution image appears on the computer's display.

## Low X-Ray Dosage

This unit's largest field of view of 170 x 120 mm (the entire head/neck area) offers dosage 25% less than an FMX, F-speed film.




Isotropic cubic voxel  
(3D Accuitomo)





Anisotropic voxel  
(Conventional CT)

## Isotropic Cubic Voxel


A voxel is the minimum unit of 3D data. 3D Accuitomo's voxel is an isotropic cube that produces images with equally fine detail in all three dimensions and minimizes artifacts produced by slice pitch and helical pitch, therefore resolution is never degraded by re-slicing. Rectangular voxels, used for conventional CT imaging (and some dental cone beam units), result in image degradation when re-sliced.



Enlarged yet smooth and distortion-free  
80  $\mu$ m high resolution image



Temporal bone:  $\phi$  60 x H 60 mm



Paranasal sinuses:  $\phi$  100 x H 100 mm

### Flat Panel Detector (FPD) Position Adjustment

Adjusting the position of the FPD reduces X-ray dosage, provides higher resolution, and minimizes distortion.

For regions such as 140 x 100 mm, moving the FPD slightly farther away from the center of the exposure area results in a better orthographic projection, which reduces distortion and improves resolution.


Optimizing collimation of the beam, depending on the size of the area, also reduces X-ray dosage and X-ray scattering as well.

### Nine Fields of View

There are nine sizes for exposure regions with diameters ranging from 40 to 170 mm.

Diameter x Height (mm)

- $\phi$  170 x H 120 mm     $\phi$  170 x H 50 mm
- $\phi$  140 x H 100 mm     $\phi$  140 x H 50 mm
- $\phi$  100 x H 100 mm     $\phi$  100 x H 50 mm
- $\phi$  80 x H 80 mm
- $\phi$  60 x H 60 mm
- $\phi$  40 x H 40 mm




Observation of the paranasal sinuses together with right and left temporal bone:  $\phi$  170 x H 50 mm

### Resolution Stays High for Even Large Areas

Resolution stays high and distortion is minimized for all regions from the smallest to the largest field of view.




Observation of the paranasal sinuses together with right and left temporal bone:  $\phi$  170 x H 120 mm



## NINE FIELDS OF VIEWS

**3D ACCUITOMO 170 OFFERS A WIDE FOV RANGE** from a couple teeth up to the entire head and neck area. By closely matching the FOV to the region of interest, patient dose is kept to a minimum.




## HIGH QUALITY 3D IMAGES WITH LOW DOSAGE

### USING A HIGH-SENSITIVITY, HIGH-RESOLUTION

flat panel detector, extremely detailed, high resolution images can be obtained of the head and neck region such as the temporomandibular joint (TMJ), erupting or impacted teeth, dental anomalies and pathology, tumors or cysts, fractures, sinuses, and cranial base for a wide range of multi-purpose diagnostic scanning.



### Spatial Resolution\* MTF: Modulation Transfer Function



*This function is based on data from a typical product.*



### Flat Panel Detector (FPD)

FPD conversion of X-ray exposure into a digital signal results in a dramatic improvement in image quality and a reduction in X-ray dosage. The FPD is not affected by magnetic fields and has superb sensitivity and resolution to produce superior 3D images with minimal distortion and a wide dynamic range expressed with a rich distribution of grayscale values.

X-rays are converted into visible light by the directly deposited CsI scintillator and then the light is converted into an electrical signal by a photo diode. The FPD is quite thin and has a long working life.

### High Resolution

Detailed images have a resolution of at least 2.0 lp/mm (MTF 10%) with a voxel size of 80  $\mu$ m.



Highly detailed imaging

### Minimal Distortion

The flatness of the detector minimizes distortion. This eliminates the need for a compensation process, as used with analog systems, to correct for distortion prior to image reconstruction.

### Wide Dynamic Range\*\*

The FPD has a wide dynamic range of 14 bit data (64 times 8 bit data). This produces a precise grayscale differentiation.

\* Spatial resolution refers to how distinct an image appears the smaller it becomes; it measures the fineness of an image. Spatial frequency is the unit of measurement of line pairs per distance (mm). As the map scale decreases, the patterns of contrast become harder to see. This is known as MTF (Modulation Transfer Function). It represents the number of line pairs per 1 mm that can be distinguished based on contrast. It is said that humans can only differentiate about 10%.

\*\* Dynamic range: Numerical values express the reproducibility of the signal and the ratio of the largest and smallest input values in dBs. The dynamic range of the digital signal is also sometimes expressed in bits. The highest signal level is taken to be the level remaining after subtracting the noise level. The value of the dynamic range indicates how weak of a signal can be reproduced, or, in other words, how high the contrast resolution will actually turn out to be.



# 3D Accu<sup>i</sup>tomo

XYZ Slice View Tomograph

READY SCAN time 1.0 sec

MODE			
CT	<b>SCOUT</b>		
SCAN			
180°	360°		
SCAN MODE			
Std	Hi-Fi	16-Res	Hi-Speed

FOV  
**Φ170 × 120**

X	0.0mm
Y	+20.0mm
Z	+122mm

**80 kV 2.0 mA**



Tube Voltage	Tube Current	Param. Memory
--------------	--------------	---------------

J. MORITA MFG. CORP.  
Kyoto Japan



IN/OUT READY BEAM

LOCK/FREE





## SIMPLE, ACCURATE POSITIONING

**THE THREE POSITIONING LASER BEAMS AND AN LCD** make patient positioning easy. The chinrest stabilizes the patient's head to avoid movement. Scout images enable even more accurate positioning.


### Easy as One, Two, Three

First, the patient's initial position is set and recorded using the three positioning laser beams. Then, the region of interest is aligned in the LCD. The chair automatically moves into the optimal position. During the X-ray exposure, the patient is stabilized by the chinrest and the headrest.


### 2-Directional Scout

For even more accurate positioning, scout images can be utilized. After positioning, two still X-ray images of coronal and sagittal views can be taken to confirm that the position is accurate. If adjustment is necessary, positioning can be changed by dragging the cursor on the monitor and moving it to the center of the region of interest. Hitting the "ready" key will then automatically move the chair, and thus the region of interest, to the desired position.


The Scout exposure (80 kV and 2.0 mA) will increase the total X-ray dosage of a Standard Mode CT exposure (90 kV and 5.0 mA) by about 2%.



3D image region of interest is well centered.



Temporal bone:  $\phi$  60 x H 60 mm. Voxel size: 125  $\mu$ m



Paranasal sinuses:  $\phi$  170 x H 120 mm. Voxel size: 250  $\mu$ m


## LIMITED CT IMAGE AREA FOR REDUCED X-RAY DOSAGE

19

**LIMIT THE X-RAY DOSAGE** by accurately determining the minimum region of interest before exposing the patient to the higher dosage CT scan.


### Imaging Area: Diameter x Height (mm)

ø 170 x H 120 mm	ø 170 x H 50 mm
ø 140 x H 100 mm	ø 140 x H 50 mm
ø 100 x H 100 mm	ø 100 x H 50 mm
ø 80 x H 80 mm	
ø 60 x H 60 mm	
ø 40 x H 40 mm	



Scout Image

## ZOOM RECONSTRUCTION WITH 80 $\mu\text{m}$ VOXEL RESOLUTION




Select a region of interest such as the temporal bone or paranasal sinus and zoom in with 80  $\mu\text{m}$  voxel resolution for a more detailed observation.

*Auditory ossicular chain zoom reconstruction*




*Ethmoid sinus zoom reconstruction*



*Auditory ossicular chain*



Volume rendering produces a detailed 3D view of internal structures.



*Lamina Cribrosa*



**RECENTLY UPDATED, i-Dixel 2.0 IMAGE PROCESSING SOFTWARE is easy to use with several added features for advanced image manipulation and observation.**



### Compatibility

i-Dixel 2.0 comes standard with 3D Accuitomo 170 and all Morita 3D units. It is compatible with Microsoft Windows® 7 and can be seamlessly integrated into a network environment.

### Version 2.0 Updates

An improved user interface makes i-Dixel 2.0 easier to use and navigate. Updated filters between i-Dixel and One Volume Viewer (viewing software) offer uniform image quality and resolution. Several new features have also been added for advanced image manipulation including real time re-slice in the highest resolution, implant planning, and a nerve canal tracing function. i-Dixel 2.0 is also fully integrated with TDO Practice Management Software – the most comprehensive and widely used endodontic software in the world. TDO users can seamlessly incorporate an Accuitomo unit into their practice with full compatibility.

### Other Key Features

- XYZ View Windows
- Re-slice
- Zoom
- Rotate
- Histogram
- Edge Enhancement
- Distance and Angle Measurement
- Negative Image
- Mirror Image
- Slice Distance Measurement
- Surface Rendering
- DICOM 3.0 Compatible
- Brightness Conversion
- Spatial Frequency Filter
- Patient Orientation Display
- Density Measurement
- Nerve Canal Tracing
- Report Comments for Images

## FEATURES

### Volume Rendering

Volume rendering of CT data produces three dimensional images. Select the area of interest and adjust the controls for the histogram to create a detailed image of very fine structures.

### Curved MPR (cMPR)

cMPR allows observation of an orthogonal representation of the dental arch or any arbitrary curve.

## SHARING IMAGE DATA

**INSTALLING i-Dixel SOFTWARE** on all intra-clinic computers enables sharing of image data on each linked client computer. Observation of images on non-network computers can be achieved with **One Data Viewer** and **One Volume Viewer** without installing i-Dixel.

### **One Data Viewer & One Volume Viewer**



These unique Morita applications allow you to view three dimensional images and volume rendered images even if the computer does not have i-Dixel software installed.

CT data can be exported from the i-Dixel application and later stored on a CD, DVD, or removable media. This information can then be transferred to a computer outside the clinic to view CT images, volume rendered images, and patient information.

Additional functions include zoom, black and white reverse, brightness and contrast adjustment, as well as optional length and angle measurement capabilities.

### **i-Dixel conforms to the following DICOM standards:**

1. Modality worklist management service class
2. Storage service class
3. Modality performed procedure step service class
4. Print management service class






*One Volume Viewer*




OSTEOMYELITIS


Osteomyelitis of the left condyle and the temporomandibular fossa.




## NASAL SINUSITIS AND INVERTED SUPERNUMERARY TOOTH

All paranasal sinuses are fully developed. Findings of the sinuses are normal except the left sinus in which there is a polyp on the frontal wall. In the left nasal cavity, the mucous membrane of concha nasalis inferior and media is thickened. The nasal septum deviates to the left. These changes are narrowing the middle meatus and the nasal air passage. The image confirms that there is no obstruction in the middle meatus or the infundibulum.







In addition, an inverted supernumerary tooth exists in the middle of the palate. It does not represent follicle enlargement or other pathological changes.



## LABIAL BONE

Labial bone graft of the upper anterior region.



## LARGE FOLLICULAR CYST


Large follicular cyst in the left body of the mandible.



## MAXILLARY SINUS ASYMMETRY

29

Marked asymmetry between right and left maxillary sinus.




## RESIDUAL ROOT

Inflammatory bone destruction associated with a residual root tip.



## IMPLANT PLANNING


Volume of the alveolar bone of the region of missing molar.



# TECHNICAL SPECIFICATIONS



## Dimensions



- 1 Outlet of computer cable and operation cable
- 2 Outlet of power supply
- 3 Outlet of operation cable
- 4 Outlet of computer cable

## Specifications

Trade Name	3D AccuTom													
	XYZ Slice View Tomograph													
Model	MCT-1													
Type	EX1/2 F17													
Input Voltage	100 / 110 / 120 VAC 220 / 230 / 240 VAC													
Power Consumption	max. 2.0 kVA													
X-ray Head														
Tube Voltage	60 - 90 kV													
Tube Current	1-10 mA (Max 8 mA : Hi-Fi, Hi-Res Mode)													
Focal Spot Size	0.5 mm													
Exposure Time (360° / 180°)	Std Mode: 17.5 / 9.0 sec. Hi-Fi Mode: 30.8 / 15.8 sec. Hi-Res Mode: 30.8 / 15.8 sec. Hi-Speed Mode: 10.5 / 5.4 sec.													
Field of View (Diameter x Height)	<table border="0"> <tr> <td>ø 170 x H 120 mm</td> <td>ø 170 x H 50 mm</td> </tr> <tr> <td>ø 140 x H 100 mm</td> <td>ø 140 x H 50 mm</td> </tr> <tr> <td>ø 100 x H 100 mm</td> <td>ø 100 x H 50 mm</td> </tr> <tr> <td>ø 80 x H 80 mm</td> <td></td> </tr> <tr> <td>ø 60 x H 60 mm</td> <td></td> </tr> <tr> <td>ø 40 x H 40 mm</td> <td></td> </tr> </table>		ø 170 x H 120 mm	ø 170 x H 50 mm	ø 140 x H 100 mm	ø 140 x H 50 mm	ø 100 x H 100 mm	ø 100 x H 50 mm	ø 80 x H 80 mm		ø 60 x H 60 mm		ø 40 x H 40 mm	
ø 170 x H 120 mm	ø 170 x H 50 mm													
ø 140 x H 100 mm	ø 140 x H 50 mm													
ø 100 x H 100 mm	ø 100 x H 50 mm													
ø 80 x H 80 mm														
ø 60 x H 60 mm														
ø 40 x H 40 mm														
Voxel Size	80 µm, 125 µm, 160 µm, 200 µm, 250 µm													
Outer Dimensions														
Main Unit (W x D x H)	1,620 mm x 1,250 mm x 2,080 mm (63-3/4" x 49-1/4" x 82")													
Control Box (W x D x H)	96 mm x 40 mm x 115 mm (3-3/4" x 1-5/8" x 4-1/2")													
Weight	Approx. 400 kg. (882 lbs)													

*X-ray protection should be provided for the patient when X-rays are emitted.  
Design and specifications are subject to change without notification.*

Images provided by:  
Fukushima Medical University  
Technical Assistance: NUBIC (Nihon University Business, Research,  
and Intellectual Property Center)  
Collaborative Development: J. Morita Corp. & Nihon University

Erkki Tammisalo  
Emeritus professor of Oral Radiology  
Tomodent, Private Laboratory of Oral Diagnostic Imaging  
Turku, Finland



## J. MORITA MFG. CORP.

Distributed by

**J. MORITA USA, INC.**

9 Mason, Irvine, CA 92618 U.S.A.  
Toll Free: (800) 831-3222, Fax: (949) 465-1095, [www.morita.com/usa](http://www.morita.com/usa)

**J. MORITA EUROPE GMBH**

Justus-von-Liebig-Strasse 27A, D-63128 Dietzenbach, Germany  
Tel: +49-6074-836-0, Fax: +49-6074-836-299, [www.morita.com/europe](http://www.morita.com/europe)

**J. MORITA AUSTRALIA & NEW ZEALAND**

247 Coward Street, Suite 2.05, Mascot NSW 2020, Australia  
Tel: +61 2 9667 3555, Fax: +61 2 9667 3577, [www.morita.com/anz](http://www.morita.com/anz)

**J. MORITA MIDDLE EAST**

Saraya Al Wessam Bldg., Tower A, Appt. 902  
Cross roads Tag Roasaa & Abu Heif, Saba Pacha, Alexandria, Egypt  
Tel: 203 58 222 94, Fax: 203 58 222 96, [www.jmoritamiddleeast.com](http://www.jmoritamiddleeast.com)

**J. MORITA CORPORATION**

33-18, 3-Chome, Tarumi-cho Suita City, Osaka, 564-8650 Japan  
Tel: +81-6-6380-2525, Fax: +81-6-6380-0585, [www.morita.com/asia](http://www.morita.com/asia)

Developed and Manufactured by

**J. MORITA MFG. CORP**

680 Higashihama Minami-cho, Fushimi-ku, Kyoto, 612-8533 Japan  
Tel: +81-75-611-2141, Fax: +81-75-622-4595, [www.morita.com/dental](http://www.morita.com/dental)